









ProVision Contact Lens Product Guide

	Product Name & Product Order Code	Product Description	Key Benefits	Power Range
Silicone Hydrogel	 ProView ELITE® Order code: PV3 (30)/PV9 (90)	DAILY MyDay® daily disposable with Smart Silicone® chemistry, stenfilcon A 54% water content, 30 Pack BC 8.4mm / DIA 14.2mm	<ul style="list-style-type: none"> • With oxygen-efficient Smart Silicone® chemistry, ProView ELITE lenses deliver the optimal balance of breathability, comfort and easy handling • Aspheric optics to minimise spherical aberration • UV inhibitor¹ 	+0.25 to +5.00D (0.25D) +5.50 to +8.00D (0.50D) -0.25 to -6.00D (0.25D) -6.50 to -12.00D (0.50D)
	 ProView ELITE® ASTIGMATISM Order code: PX1 (30)/PX2 (90)	DAILY MyDay® daily disposable toric with Smart Silicone® chemistry & Optimised Toric Lens Geometry™, stenfilcon A 54% water content, 30 Pack BC 8.6mm / DIA 14.5mm	<ul style="list-style-type: none"> • With oxygen-efficient Smart Silicone® chemistry, ProView ELITE ASTIGMATISM lenses deliver the optimal balance of breathability, comfort and easy handling • Optimised Toric Lens Geometry™ for excellent lens stability; wide ballast band, uniform horizontal ISO thickness, and smooth continuous surface • UV inhibitor¹ 	+0.50 to +6.00DS (0.50D)* Plano to -6.00DS (0.25D) -6.50 to -10.00DS (0.50D) Cyl -0.75, -1.25, -1.75, -2.25DC 10, 20, 70, 80, 90, 100, 110, 160, 170, 180 for -0.75 to -1.75DC 10, 20, 90, 160, 170, 180 for -2.25DC
	 ProView PERFORM® Order code: 8JE (30)/8JH (90)	DAILY clariti® 1 day, somofilcon A 56% water content, 30 Pack BC 8.6mm / DIA 14.1mm	<ul style="list-style-type: none"> • Silicone hydrogel at the price of a hydrogel • A highly wetttable lens surface, supporting excellent all-day comfort • Aspheric optics to minimise spherical aberration • UV inhibitor¹ 	+0.50 to +6.00D (0.25D) +6.50 to +8.00D (0.50D) -0.50 to -6.00D (0.25D) -6.50 to -10.00D (0.50D)
	 ProView PERFORM® ASTIGMATISM Order code: 8JF (30)/8Ji (90)	DAILY clariti® 1 day toric, somofilcon A 56% water content, 30 Pack BC 8.6mm / DIA 14.3mm	<ul style="list-style-type: none"> • Silicone hydrogel at the price of a hydrogel • Aspheric front optic zone to minimise spherical aberration • UV inhibitor¹ • Easy fitting with OptiExpert® app 	Plano to +4.00DS (0.25D) -0.25 to -6.00DS (0.25D) -6.50 to -9.00DS (0.50D) Cyl -0.75, -1.25, -1.75, -2.25DC 10, 20, 60, 70, 80, 90, 100, 110, 120, 160, 170, 180 for -0.75 to -1.75DC (all minus) 20, 70, 90, 110, 160, 180 for -0.75 to -1.75DC (all plus) , 10, 20, 90, 160, 170, 180 for -2.25DC (all minus)
	 ProView PERFORM® MULTIFOCAL Order code: 8JG (30)/8Ji (90)	DAILY clariti® 1 day multifocal, somofilcon A 56% water content, 30 Pack BC 8.6mm / DIA 14.1mm	<ul style="list-style-type: none"> • Silicone hydrogel at the price of a hydrogel • Unique dual-zone intermediate optics with low and high add powers • UV inhibitor¹ • Easy fitting with OptiExpert® app 	Plano to +5.00D (0.25D) -0.25 to -6.00D (0.25D) Low: for spec adds up to +2.25D High: for spec adds over +2.25D
Hydrogel	 ProView® Order code: ZPP (30)/ZP9 (90)	DAILY Proclear® 1 day, omafilcon A 60% water content, 30 Pack BC 8.7mm / DIA 14.2mm	<ul style="list-style-type: none"> • Patented PC Technology™ - allows the lens to stay hydrated, remaining fresh and comfortable all day • Aspheric optics to reduce spherical aberration 	+0.25 to +5.00D (0.25D) +5.50 to +8.00D (0.50D) -0.25 to -6.00D (0.25D) -6.50 to -12.00D (0.50D)
	 ProView® MULTIFOCAL Order code: PVF	DAILY Proclear® 1 day multifocal, omafilcon A 60% water content, 30 Pack BC 8.7mm / DIA 14.2mm	<ul style="list-style-type: none"> • Patented PC Technology™ - allows the lens to stay hydrated, remaining fresh and comfortable all day • Aspheric optics to reduce spherical aberration • Easy fitting with OptiExpert® app 	Plano to +6.00D (0.25D) -0.25 to -6.00D (0.25D) -6.50 to -10.00D (0.50D)

Product Name	Product Description	Key Benefits	Power Range
	ProActive AIR® Order code: 6PO FORTNIGHTLY cAir®, fanfilcon A 55% water content, 6 Pack BC 8.4mm / DIA 14.2mm	<ul style="list-style-type: none"> The latest 2-weekly silicone hydrogel technology High oxygen transmissibility and water content UV Inhibitor¹ 	+0.25 to +6.00D (0.25D) +6.50 to +8.00D (0.50D) -0.25 to -6.00D (0.25D) -6.50 to -12.00D (0.50D)
	ProActive AIR® ASTIGMATISM Order code: PROAVT FORTNIGHTLY cAir® toric, fanfilcon A 55% water content, 6 Pack BC 8.5mm / DIA 14.5mm	<ul style="list-style-type: none"> The latest 2-weekly silicone hydrogel technology Optimised Toric Lens Geometry™ for excellent lens stability; wide ballast band, uniform horizontal ISO thickness, and smooth continuous surface UV Inhibitor¹ Replacement for ProActive ASTIGMATISM² 	Plano to +6.00DS (0.25D) +6.50 to +8.00DS (0.50D) -0.25 to -6.00DS (0.25D) -6.50 to -10.00DS (0.50D) Cyl -0.75, -1.25, -1.75, -2.25DC Axis 10 to 180 (10° steps)
	ProActive® ASTIGMATISM Order code: ZPT FORTNIGHTLY Avaira® toric, enfilcon A 46% water content, 6 Pack BC 8.5mm / DIA 14.5mm	<ul style="list-style-type: none"> Aquaform® Technology- unique material technology creates a naturally wettable silicone hydrogel contact lens Optimised Toric Lens Geometry™ for excellent lens stability; wide ballast band, uniform horizontal ISO thickness, and smooth continuous surface UV Inhibitor¹ 	Plano to +6.00DS (0.25D) -0.25 to -6.00DS (0.25D) -6.50 to -10.00DS (0.50D) Cyl -0.75, -1.25, -1.75, -2.25DC Axis 10 to 180 (10° steps)
	ProBalance® ENHANCE Order code: APC MONTHLY Biofinity Energys™, comfilcon A 48% water content, 3 Pack BC 8.6mm / DIA 14.0mm	<ul style="list-style-type: none"> Digital Zone Optics™ lens design - helps to ease accommodative burden without impacting distance visual acuity Aquaform® Technology - unique material technology creates a naturally wettable silicone hydrogel contact lens Daily, flexible or up to 29 nights' continuous wear³ 	Plano to +6.00D (0.25D) +6.50 to +8.00D (0.50D) -0.25 to -6.00D (0.25D) -6.50 to -12.00D (0.50D)
	ProBalance® Order code: APB MONTHLY Biofinity®, comfilcon A 48% water content, 3 Pack BC 8.6mm / DIA 14.0mm	<ul style="list-style-type: none"> Aquaform® Technology - unique material technology creates a naturally wettable silicone hydrogel contact lens. Soft, flexible low modulus lens with high levels of oxygen transmissibility Daily, flexible or up to 29 nights' continuous wear³ 	+0.25 to +6.00D (0.25D) +6.50 to +8.00D (0.50D) -0.25 to -6.00D (0.25D) -6.50 to -12.00D (0.50D)
	ProBalance® ASTIGMATISM Order code: APG MONTHLY Biofinity® toric, comfilcon A 48% water content, 3 Pack BC 8.7mm / DIA 14.5mm	<ul style="list-style-type: none"> Optimised Toric Lens Geometry™ for excellent lens stability; wide ballast band, uniform horizontal ISO thickness, and smooth continuous surface Aquaform® Technology - unique material technology creates a naturally wettable silicone hydrogel contact lens Easy fitting with OptiExpert® app 	Plano to +6.00DS (0.25D) +6.50 to +8.00DS (0.50D) -0.25 to -6.00DS (0.25D) -6.50 to -10.00DS (0.50D) Cyl -0.75, -1.25, -1.75, -2.25DC Axis 10 to 180 (10° steps)
	ProBalance® MULTIFOCAL Order code: ZPD (D)/APN (N) MONTHLY Biofinity® multifocal, comfilcon A 48% water content, 3 Pack BC 8.6mm / DIA 14.0mm	<ul style="list-style-type: none"> Aquaform® Technology - unique material technology creates a naturally wettable silicone hydrogel contact lens Balanced Progressive™ Technology delivers excellent vision at any distance A streamlined fitting approach and four add powers Centre-distance and centre-near designs Easy fitting with OptiExpert® app 	Plano to +6.00D (0.25D) -0.25 to -6.00D (0.25D) -6.50 to -10.00D (0.50D) Design: D and N lens Add power +1.00, +1.50, +2.00, +2.50D

Silicone Hydrogel

1. Warning: UV-absorbing contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, as they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.
 2. ProActive ASTIGMATISM is to be replaced by ProActive AIR ASTIGMATISM 3. Biofinity lenses are recommended for monthly replacement and are CE marked under the Medical Devices Directive for occasional overnight or up to 29 nights' continuous wear under the guidance of the eyecare practitioner.
 * Plus powers available 2018